

Highcon PCR Purification Kits



Contents

Kit Contents	2
Principle	3
Important Notes	3
Highcon PCR Purification Kit Protocol	4
Using a Microcentrifuge	4
Using a Vacuum Manifold	5
Troubleshooting Guide	6
Ordering Information	7

Kit Contents

Catalog No.	2315	2312
Number of preparations	50	200
Nucleogen Highcon Spin Columns	50	200
Collection tubes (2 ml)	50	200
PCR Purification Buffer*	+	+
Wash A' Solution (concentrate)	+	+
Elution Buffer	+	+

* PCR Purification Buffer contain chaotropic salts which are irritants and not compatible with disinfecting agents containing bleach. Take appropriate laboratory safety measures and wear gloves when handling.

Principle

The Nucleogen PCR Purification Buffer is optimized for efficient recovery of DNA and removal of contaminants. DNA adsorbs to the silica-membrane in the presence of high salt while contaminants pass through the column. Impurities are efficiently washed away, and the pure DNA is eluted with Elution Buffer or Water. The Highcon PCR Purification columns offer two handling options — as an alternative to processing the Highcon columns in a microcentrifuge, they can now also be used on any commercial vacuum manifold with luer connectors.

Adsorption to Highcon PCR Purification membrane

The Highcon PCR Purification membrane is uniquely adapted to isolate DNA from aqueous solutions, and up to 7 ug DNA can bind to Highcon PCR Purification column.

Washing

During the DNA adsorption step, unwanted primers and impurities, such as salts, enzymes, unincorporated nucleotides, agarose, dyes, ethidium bromide, oils, and detergents (e.g., DMSO, Tween® 20) do not bind to the membrane, but flow through the column. Salts are quantitatively washed away by the ethanol-containing Wash A' Solution. Any residual Wash A' Solution, which may interfere with subsequent enzymatic reactions, is removed by an additional centrifugation step.

Elution in low-salt solutions

Elution efficiency is strongly dependent on the salt concentration and pH of the Elution Buffer. Contrary to adsorption, elution is most efficient under basic conditions and low salt concentrations. DNA is eluted with 9 or 10 ul of the provided Elution Buffer (10 mM Tris-Cl, pH 8.5), or Water. The maximum elution efficiency is achieved between pH 7.0 and 8.5.

When using water to elute, make sure that the pH is within this range. In addition, DNA must be stored at -20°C when eluted with water since DNA may degrade in the absence of a buffering agent. Elution with TE (10 mM Tris-Cl, 1 mM EDTA, pH 8.0) is possible, but not recommended because EDTA may inhibit subsequent enzymatic reactions.

DNA yield and concentration

DNA yield depends on the following three factors: the volume of Elution Buffer, how the buffer is applied to the Highcon column, and the incubation time of the buffer on the Highcon column. 20–30 ul of elution buffer completely covers the Highcon PCR Purification membrane, ensuring maximum yield, even when not applied directly to the center of the membrane. Elution with 9-10 ul requires the buffer to be added directly to the center of the membrane.

Important Notes

Please read the following notes before starting any of the Nucleogen Highcon PCR Purification procedures.

Before equipment

- **Add ethanol(100%) to Wash A' Solution** before use (see bottle label for volume).
- All centrifuge steps are at 13,000rpm (~17,900 x g) in a conventional tabletop microcentrifuge.

Vacuum notes:

- Switch off vacuum between steps to ensure that a consistent, even vacuum is applied during manipulations.
- Wear safety glasses when working near a manifold under pressure.
- The vacuum pressure is the pressure differential between the inside of the manifold and the atmosphere

Highcon PCR Purification Kit Protocol Using a Microcentrifuge

Please read “Important Notes” on pages 3 before starting.

This protocol is designed to purify single- or double-stranded DNA fragments from PCR and other enzymatic reactions. For cleanup of other enzymatic reactions. Fragments ranging from 70 bp to 10 kb are purified from primers, nucleotides, polymerases, and salts using Highcon spin columns in a microcentrifuge.

- 1. Add 5 volumes of PCR Purification Buffer to 1 volume of the PCR sample and mix.**
For example, add 500 ul of PCR Purification Buffer to 100 ul PCR sample (not including oil).
- 2. Place a Highcon column in a provided 2 ml collection tube.**
- 3. To bind DNA, apply the sample to the Highcon column and centrifuge for 30-60 sec.**
- 4. Discard flow-through and place the Highcon column back in the same tube.**
Collection tubes are re-used to reduce plastic waste.
- 5. To wash, add 750 ul Wash A' Solution to the Highcon column and centrifuge for 30-60 sec.**
- 6. Discard flow-through and place the Highcon column back in the same tube.**
Centrifuge the Highcon column for an additional 1 min.
IMPORTANT: Residual ethanol from Wash A' Solution will not be completely removed unless the flow-through is discarded before this additional centrifugation.
- 7. Place Highcon column in a clean 1.5 ml microcentrifuge tube.**
- 8. To elute DNA, add 10 ul Elution Buffer (10 mM Tris-Cl, pH 8.5) or H₂O to the center of the Highcon column membrane, let the Highcon column stand for 1 min and then centrifuge for 1 min.**
IMPORTANT: Ensure that the elution buffer is dispensed directly onto the center of the membrane for complete elution of bound DNA. The average eluate volume is 9 ul from 10 ul elution buffer volume.
Elution efficiency is dependent on pH. The maximum elution efficiency is achieved between pH 7.0 and 8.5.

Using a Vacuum Manifold

Please read “Important Notes” on pages 3 before starting.

- 1. Add 5 volumes of PCR Purification Buffer to 1 volume of the PCR sample and mix.**
For example, add 500 ul of PCR Purification Buffer to 100 ul PCR sample (not including oil).
- 2. Prepare the vacuum manifold and Highcon PCR Purification columns:**
 - Insert each Highcon column into a luer connector on the Luer Adapter(s) in the manifold.
- 3. To bind DNA, load the samples into the Highcon columns by decanting or pipetting, and apply vacuum. After the samples have passed through the Highcon column, switch off the vacuum source.**
The maximum loading volume of the highcon column is 800 ul. For sample volumes greater than 800 ul simply load again.
- 4. To wash, add 750 ul of Wash A' Solution to each Highcon column and apply vacuum.**
- 5. Transfer each Highcon column to a microcentrifuge tube or the provided 2 ml collection tubes.**
Centrifuge for 1 min at 13,000 rpm (~17,900 x g).
IMPORTANT: This spin is necessary to remove residual ethanol (Wash A' solution).
- 6. Place each Highcon column into a clean 1.5 ml microcentrifuge tube.**
- 7. To elute DNA, add 10 ul Elution Buffer (10 mM Tris-Cl, pH 8.5) or H₂O to the center of the Highcon column membrane, let the Highcon column stand for 1 min and then centrifuge for 1 min.**
IMPORTANT: Ensure that the elution buffer is dispensed directly onto the center of the membrane for complete elution of bound DNA. The average eluate volume is 9 ul from 10 ul elution buffer volume.
Elution efficiency is dependent on pH. The maximum elution efficiency is achieved between pH 7.0 and 8.5.

Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise.

Comments and suggestions

Low or no yield

- | | |
|---|--|
| a) Wash A' Solution did not contain ethanol | Ethanol must be added to Wash A' Solution (concentrate) before use. Repeat procedure with correctly prepared Wash A' Solution. |
| b) Inappropriate elution buffer | DNA will only be eluted efficiently in the presence of low-salt buffer (e.g., Elution Buffer: 10 mM Tris-Cl, pH 8.5) or water. |
| c) Elution Buffer
Incorrectly dispensed | Add Elution Buffer to the center of the membrane to ensure that the buffer completely covers the membrane.
This is particularly important when using small elution volumes (10 ul). |

DNA does not perform well, e.g., in ligation reactions

- | | |
|--|---|
| a) Salt concentration in eluate too high | Modify the wash step by incubating the column for 5 min at room temperature after adding 750 ul of Wash A' Solution, then centrifuge. |
| b) Elution contains residual ethanol | Ensure that the wash flow-through is drained from the collection tube and that the Highcon column is then centrifuged at 13,000 rpm (~17,900 x g) for an additional 1 min. |
| c) Elution contains primer-dimers | Primer-dimers formed are longer than 20 bp, and are not completely removed. After the binding step, wash the Highcon column with 750 ul of a 35% guanidine hydrochloride aqueous solution (35 g in 100 ml). Continue with the Wash A' Solution wash step and the elution step as in the protocol. |
| d) Eluate contains denatured ssDNA, which appears as smaller smeared band on an analytical gel | Use the eluted DNA to prepare the subsequent enzymatic reaction but omit the enzyme. To reanneal the ssDNA, incubate the reaction mixture at 95°C for 2 min, and allow the tube to cool slowly to room temperature. Add the enzyme and proceed as usual. Alternatively, the DNA can be eluted in 10 mM Tris buffer containing 10 mM NaCl. The salt and buffering agent promote the renaturation of DNA strands. However the salt concentration of the eluate must then be considered for subsequent applications. |

Ordering Information

Products	Contents	Cat. No.	
RNA Stabilization Reagent	Tube (50 x 1.5 ml)	3502	
	Tube (20 x 5 ml)	3205	
	100 ml	3100	
	250 ml	3250	
Plasmid Purification Mini Kit	(200) for negative strain	200 preps	5112
	(200) for positive strain	200 preps	7112
Plasmid Purification Midi Kit	(10)	10 preps	6101
	(50)	50 preps	6105
	(100)	100 preps	6110
Plasmid Purification Maxi Kit	(6)	6 preps	7106
	(24)	24 preps	7124
	(50)	50 preps	7150
Gel Extraction Kit	(50)	50 preps	5215
	(200)	200 preps	5212
Highcon Gel Extraction Kit	(50)	50 preps	2215
	(200)	200 preps	2212
Bead Type (Nal) Gel Extraction Kit	(200)	200 preps	1232
	(400)	400 preps	1234
	(600)	600 preps	1236
PCR Purification Kits	(50)	50 preps	5315
	(200)	200 preps	5312
Highcon PCR Purification Kit	(50)	50 preps	2315
	(200)	200 preps	2312
DNA Clean-up Kits	(50)	50 preps	1415
	(200)	200 preps	1412
Genomic Blood Spin Mini Kit	(50)	50 preps	1515
	(200)	200 preps	1512
Genomic Blood Spin Midi Kit	(20)	20 preps	6520
	(50)	50 preps	6550
	(100)	100 preps	6500
Genomic Blood Spin Maxi Kit	(6)	6 preps	7506
	(24)	24 preps	7524
	(50)	50 preps	7550

Ordering Information

Products		Contents	Cat. No.
Genomic Cell / Tissue Spin Mini Kit	(50)	50 preps	1545
	(200)	200 preps	1542
Genomic Cell / Tissue Spin Midi Kit	(20)	20 preps	
	(50)	50 preps	
	(100)	100 preps	
Genomic DNA Isolation, Flexible		100 Isolation	1521
		500 Isolation	1525
		10 ml x 100 Isolation	
Apoptotic DNA Ladder Kit		50 preps	2505
96 PCR Purification Kit			
	4 x 96 plates(binding, elution), buffer, tape		4304
	25 x 96 plates(binding, elution), buffer, tape		4325
	50 x 96 plates(binding, elution), buffer, tape		2 x 4325
96 Plasmid Purification Kit			
	4 x 96 plates(clarification, binding, elution), buffer, tape		4104
	25 x 96 plates(clarification, binding, elution), buffer, tape		4125
96 Genomic Blood Spin Kit			
	4 x 96 plates(binding, elution), buffer, tape		
	25 x 96 plates(binding, elution), buffer, tape		
96 Genomic Cell / Tissue Spin Kit			
	4 x 96 plates(binding, elution), buffer, tape		
	25 x 96 plates(binding, elution), buffer, tape		
	50 x 96 plates(binding, elution), buffer, tape		2 x



www.nucleogen.com

Tel) 080-484-2007

031-315-6644

Fax) 031-312-2008